



# PDO TENSIO

Designed for professionals

uarlino TENSIO.indd 1 06/10/20 18:0

### Great Lifting Effect

## TENSIO

Threads molding manufacture, tightening force Produces a great liftting effect

Type W cannula, which allows more easily anchor in the subcutaneous tissue



#### **INDICATIONS**

Patients with medium sagging and heavy tissues

Body treatments (buttocks, breasts, abdomen)

Jaw line with great flaccidity



#### **PRESENTATION**


2 blisters per pouch with 2 threads = 4 threads (thread + cannula)
Boxes of 10 pouches, boxes of 3 pouches


CODICE	Ø AGO	LUNGH. AGO	LUNGH. FILO	USP FILO	Ø FILO
VTE181	18G	100 mm	165 mm	3	0,60 mm



#### **ESQUEMAS**







#### **KEY POINTS**

W-cannula allows safe and easy insertion into tissue, polished-treated special to generate minimal trauma.

PDO threads raw from the number one polydioxanone manufacturer, plus the thread length is 1 cm longer than other brands.

Dry manufacturing process, vacuum packed, with nitrogen gas for ensure high thread quality.

The spatial geometry of the thread allows greater glide through the cannula.

#### **BEFORE AND AFTER**









Foto courtesy of Dr.ssa Vicenta Llorca





Designed for professionals

**( €** 2292 ISO 13485



vliftpro.org

ver. 10/2020